

Zinc undecylenate


Catalog No: #S0744

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Zinc undecylenate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	431.92
Formulation	C ₂₂ H ₃₈ O ₄ Zn
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: [Zn++].OC(=O)CCCCCCCC=C
Pathways: Others
Receptor: Others
Boiling pt: 300.8°C
Melting pt: 118-121°C
Solubility: DMSO

References

1. Lin MC, et al. J Chromatogr A. 2006 Jun 30;1119(1-2):264-9.

Note: This product is for in vitro research use only