

Disodium monofluorophosphate

Catalog No: #S0745


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Disodium monofluorophosphate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Phosphatase inhibitor; Phosphorylase inhibitor; PKM inhibitor
Calculated MW	143.95
Formulation	FNa2O3P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: [O-]P([O-])(F)=O Pathways: Metabolism; Others Receptor: phosphorylase phosphatase; alkaline phosphatase; pyruvate kinase Boiling pt: Melting pt: 625°C Solubility: Water: 29 mg/mL (201.45 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Pearce EI. Caries Res, 1983, Suppl 1, 21-35.

Note: This product is for in vitro research use only