

Furazolidone

Catalog No: #S0751


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Furazolidone
Brief Description	Inhibitors
Purification	99.98%
Target Name	MAO inhibitor
Calculated MW	225.16
Formulation	C8H7N3O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: C1COC(=O)N1/N=C/C2=CC=C(O2)[N+](=O)[O-]
Pathways: Metabolism
Receptor: MAO
Boiling pt: 353.4°C
Melting pt: 254-256°C
Solubility: DMSO: > 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

1. Timperio AM, et al. Xenobiotica. 2003 Feb;33(2):153-67.

Note: This product is for in vitro research use only