

Salicylanilide


Catalog No: #S0753

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Salicylanilide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Reverse Transcriptase inhibitor
Calculated MW	213.24
Formulation	C13H11NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:OC1=CC=CC=C1C(=O)NC1=CC=CC=C1Pathways:Microbiology&virologyReceptor:Reverse TranscriptaseBoiling pt:294.3CMelting pt:136-138CSolubility:Ethanol: 18 mg/mL (84.41 mM); DMSO: 43 mg/mL (201.66 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Krtk? M, Vinsov J. Curr Pharm Des. 2011;17(32):3494-505.

Note: This product is for in vitro research use only