

Piroxicam

Catalog No: #S0758


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Piroxicam
Brief Description	Inhibitors
Purification	99.59%
Target Name	COX inhibitor
Calculated MW	331.35
Formulation	C15H13N3O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CN1C(C(=O)NC2=CC=CC=N2)=C(O)C2=CC=CC=C2S1(=O)=O
 Pathways: Immunology/Inflammation;
 Neuroscience
 Receptor: COX
 Boling pt: Melting pt: 198-200°C
 Solubility: Ethanol: 1 mg/mL (3.01 mM); DMSO: 66 mg/mL (199.18 mM)
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Blanco FJ, et al. J Rheumatol. 1999 Jun;26(6):1366-73.

Note: This product is for in vitro research use only