

Cholesterol

Catalog No: #S0760


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cholesterol
Brief Description	Inhibitors
Purification	95.00%
Target Name	MRP inhibitor; ROR
Calculated MW	386.65
Formulation	C27H46O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

systemSMILES: C[C@H](CCCC(C)C)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC=C4[C@@]3(CC[C@@H](C4)O)C)C

Pathways: Immunology/Inflammation; Metabolism Receptor: MRP1; RORα

Boling pt: 360°C Melting pt: 148-150°C Solubility: DMSO: 0.03 mg/mL (0.07 mM) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wu Y, et al. Biochem Biophys Res Commun. 2015 Feb 20;457(4):614-20.

Note: This product is for in vitro research use only