

Trans-Anethole

Catalog No: #S0762


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trans-Anethole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	148.2
Formulation	C10H12O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
SMILES: C/C=C/C1=CC=C(C=C1)OC
Pathways: Others
Receptor: Others
Boiling pt: 234-237°C
Melting pt: 20-21°C
Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Ryu J, et al. J Agric Food Chem. 2005 Jul 27;53(15):5954-8.

Note: This product is for in vitro research use only