

Diolane


Catalog No: #S0763

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Diolane
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	118.17
Formulation	C6H14O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC(O)CC(C)(C)OPathways:OthersReceptor:OthersBoiling pt:197CMelting pt:-40CSolubility:Soluble in EthanolAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Magneron I, et al. Environ Sci Technol. 2003 Sep 15;37(18):4170-81.

Note: This product is for in vitro research use only