

4-Methylpyrazole


Catalog No: #S0765

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Methylpyrazole
Brief Description	Inhibitors
Purification	99.50%
Target Name	AChR inhibitor; CAT inhibitor
Calculated MW	82.1
Formulation	C4H6N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: CC1=CN=C1
Pathways: Neuroscience; oxidation-reduction
Receptor: ADH; CAT
Boiling pt: 99-100°C
Melting pt: 13°C
Solubility: Water: 146 mg/mL
Appearance: Yellow Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Nagata H, et al. Am J Physiol. 1996 Dec;271(6 Pt 1):G1028-33.

Note: This product is for in vitro research use only