

Sulfapyridine

Catalog No: #S0767


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulfapyridine
Brief Description	Inhibitors
Purification	99.85%
Target Name	Autophagy
Calculated MW	249.29
Formulation	C ₁₁ H ₁₁ N ₃ O ₂ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: NC1=CC=C(C=C1)S(=O)(=O)NC1=CC=CC=N1
Pathways: Autophagy
Receptor: DHPS
Boiling pt: Melting pt: 191-193°C
Solubility: DMSO: 1 mg/mL (4.01 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. McDonald M, et al. J Chromatogr A. 2009 Nov 13;1216(46):8110-6.

Note: This product is for in vitro research use only