

Aceglutamide

Catalog No: #S0769


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aceglutamide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Amino Acids and Derivatives
Calculated MW	188.18
Formulation	C7H12N2O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:CC(=O)N[C@@H](CCC(=O)N)C(=O)OPathways:OthersReceptor:N-acylamino acid racemaseBoiling pt:604.9CMelting pt:206-208CSolubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Magnusson I, et al. Metabolism. 1989 Aug;38(8 Suppl 1):82-8.

Note: This product is for in vitro research use only