

Coumarin


Catalog No: #S0775

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Coumarin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Factor Xa inhibitor; Vitamin antagonist
Calculated MW	146.14
Formulation	C ₉ H ₆ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: O=C1OC2=CC=CC=C2C=C1 Pathways: Metabolism Receptor: Factor Xa; VKBoling pt: 298CMelting pt: 68-73CSolubility: Ethanol: 29 mg/mL (198.43 mM); DMSO: 29 mg/mL (198.43 mM) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Egan D, et al. Drug Metab Rev, 1990, 22(5), 503-529.

Note: This product is for in vitro research use only