

CYP2A13 Polyclonal Antibody

Catalog No: #40812

Package Size: #40812-1 50ul #40812-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CYP2A13 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	CYP2A13 Polyclonal Antibody detects endogenous levels of CYP2A13 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CYP2A13.
Target Name	CYP2A13
Other Names	CYP2A13; Cytochrome P450 2A13; CYP11A13
Accession No.	Swiss-Prot: Q16696NCBI Gene ID: 1553
Uniprot	Q16696
GeneID	1553;
SDS-PAGE MW	56kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of A549 cells using CYP2A13 Polyclonal Antibody



Note: This product is for in vitro research use only