

Phenacetin

Catalog No: #S0778


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phenacetin
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor
Calculated MW	179.22
Formulation	C10H13NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CCOC1=CC=C(NC(C)=O)C=C1Pathways:Immunology/Inflammation; NeuroscienceReceptor:COX-1Boiling pt:132CMelting pt:133-136CSolubility:Ethanol: 36 mg/mL (200.87 mM); DMSO: 36 mg/mL (200.87 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Brix AE. Toxicol Pathol. 2002 Nov-Dec;30(6):672-4.

Note: This product is for in vitro research use only