

Dithiopropanol


Catalog No: #S0781

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dithiopropanol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	124.22
Formulation	C3H8OS2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: OCC(S)CS
Pathways: Others
Receptor: Cadmium; Mercury; Others
Boiling pt: 223.4°C
Melting pt: -11.5°C
Solubility: DMSO: 10 mM
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Flora SJ, et al. Int J Environ Res Public Health. 2010, 7(7):2745-88.

Note: This product is for in vitro research use only