

Clidinium bromide

Catalog No: #S0786


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Clidinium bromide
Brief Description	Inhibitors
Purification	99.86%
Target Name	AChR antagonist
Calculated MW	432.35
Formulation	C ₂₂ H ₂₆ BrNO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism

system SMILES: [Br-].[N+]12CCCC(CC1)C(C2)OC(=O)C(O)(C1=CC=CC=C1)C1=CC=CC=C1 Pathways: Neuroscience Receptor: mAChR Boling

pt: Melting pt: 240-241°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Pathak A, et al. J Chromatogr Sci. 2010 Mar;48(3):235-9.

Note: This product is for in vitro research use only