

Butylated hydroxyanisole


Catalog No: #S0787

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Butylated hydroxyanisole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	180.25
Formulation	C ₁₁ H ₁₆ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: COC1=CC(=C(O)C=C1)C(C)(C)C
Pathways: Others
Receptor: Others
Boiling pt: 264-270°C
Melting pt: 58-60°C
Solubility: DMSO: 10 mM
Appearance: White Sheet Solid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lam LK, et al. J Med Chem. 1979 May; 22(5):569-71.

Note: This product is for in vitro research use only