

Carbaryl


Catalog No: #S0790

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carbaryl
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR inhibitor
Calculated MW	201.22
Formulation	C ₁₂ H ₁₁ NO ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CNC(=O)OC1=CC=CC2=CC=CC=C12 Pathways: Neuroscience Receptor: AChR Boiling pt: 315°C Melting pt: 142-146°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wang HP, et al. Environ Toxicol. 2014 Oct; 29(10):1193-200.

Note: This product is for in vitro research use only