

Cimetidine

Catalog No: #S0797


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Cimetidine |
| Brief Description | Inhibitors |
| Purification | 98.85% |
| Target Name | Histamine Receptor antagonist |
| Calculated MW | 252.34 |
| Formulation | C10H16N6S |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Metabolism system
 SMILES: CN=C(NCCSCC1=C(C)NC=N1)NC#N
 Pathways: GPCR/G Protein; Immunology/Inflammation;
 Neuroscience
 Receptor: H2 receptor
 Boiling pt: 488C
 Melting pt: 139-144C
 Solubility: DMSO: 51mg/ml(202.11 mM); Water:

References

1. Hernandez-Munoz R, et al. Lab Invest. 2000 Aug; 80(8): 1161-9.

Note: This product is for in vitro research use only