

## Acetrizoiic acid

Catalog No: #S0803


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Acetrizoiic acid
Brief Description	Inhibitors
Purification	96.00%
Target Name	Others
Calculated MW	556.86
Formulation	C9H6I3NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CC(=O)NC1=C(I)C=C(I)C(C(O)=O)=C1IPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. MCCONNELL F, MERSEREAU WA. J Can Assoc Radiol. 1964 Mar;15:14-22.

Note: This product is for in vitro research use only