

## Basibrol

Catalog No: #S0807


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Basibrol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	282.76
Formulation	C <sub>2</sub> H <sub>3</sub> Br <sub>3</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
SMILES: OCC(Br)(Br)Br  
Pathways: Others  
Receptor: Others  
Boiling pt: 92-93°C  
Melting pt: 73-79°C  
Solubility: Soluble in Water  
Appearance: White Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. McDowell A, et al. J Pharm Bioallied Sci. 2014 Jan;6(1):16-21.

Note: This product is for in vitro research use only