

Ethambutol dihydrochloride

Catalog No: #S0810


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ethambutol dihydrochloride
Brief Description	Inhibitors
Purification	90.00%
Target Name	Antibacterial activator
Calculated MW	277.24
Formulation	C ₁₀ H ₂₆ Cl ₂ N ₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC[C@@H](CO)NCCN[C@@H](CC)CO.Cl.Cl Pathways: Microbiology & virology Receptor: arabinosyl transferase
 Boling pt: 345.3°C Melting pt: 198-200°C Solubility: Water: 56 mg/mL (201.99 mM); DMSO: 56 mg/mL (201.99 mM) Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Tuberculosis (Edinb), 2008, 88(2), 102-105.

Note: This product is for in vitro research use only