

Spiramycin

Catalog No: #S0819


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Spiramycin
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	843.07
Formulation	C43H74N2O14
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C[C@@H]1C/C=C/C=C/[C@@H]([C@@H](C[C@@H]([C@@H]([C@H]([C@@H](CC(=O)O1)O)OC)OC2C(C(C(O2)C)OC3CC(C(C(O3)C)O)(C)O)N(C)C)O)CC(=O)C)OC4CCC(C(O4)C)N(C)C Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 100 mg/mL (118.61 mM); DMSO: 100 mg/mL (118.61 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Poulsen SM, et al. J Mol Biol, 2000, 304(3), 471-481.

Note: This product is for in vitro research use only