

Triamcinolone acetonide

Catalog No: #S0820


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Triamcinolone acetonide
Brief Description	Inhibitors
Purification	99.61%
Target Name	GR inducer
Calculated MW	434.5
Formulation	C24H31FO6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory

systemSMILES:C[C@]12C[C@@H]([C@]3([C@H]([C@@H]1C[C@@H]4[C@]2(OC(O4)(C)C)C(=O)CO)CCC5=CC(=O)C=C[C@@]53C)F)O

Pathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:Melting pt:274-278CSolubility:Ethanol: 13 mg/mL (29.91 mM); DMSO: 87 mg/mL (200.23 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nehme A, et al. BMC Med Genomics. 2009, 2:58.

Note: This product is for in vitro research use only