

Iproniazid Phosphate


Catalog No: #S0821

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Iproniazid Phosphate
Brief Description	Inhibitors
Purification	98.00%
Target Name	MAO inhibitor
Calculated MW	277.21
Formulation	C ₉ H ₁₃ N ₃ O ₄ P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: OP(O)(O)=O.CC(C)NNC(=O)C1=CC=NC=C1
 Pathways: Metabolism
 Receptor: MAO
 Boiling pt: 265.9°C
 Melting pt: 181°C
 Solubility: DMSO: 35mg/mL
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Fagervall I, et al. Biochemical pharmacology, 1986, 35 (8): 1381C7.

Note: This product is for in vitro research use only