

## Ebselen

Catalog No: #S0825


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ebselen
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor
Calculated MW	274.18
Formulation	C13H9NOSe
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Nervous systemSMILES:O=C1N([Se]C2=CC=CC=C2)C1=CC=CC=C1Pathways:Immunology/Inflammation;  
NeuroscienceReceptor:COXBoling pt:Melting pt:178-181CSolubility:Water: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility ,  
please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Schewe T. Gen Pharmacol, 2009, 26 (6): 1153C69.

Note: This product is for in vitro research use only