

Cefsulodine sodium

Catalog No: #S0827


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefsulodine sodium
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	554.52
Formulation	C ₂₂ H ₁₉ N ₄ NaO ₈ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1C(=C(N2[C@H](S1)[C@@H](C2=O)NC(=O)[C@@H](C3=CC=CC=C3)S(=O)(=O)[O-])C(=O)[O-])C[N+](=O)C=C(C=C4)C(=O)N.[Na+] Pathways: Others Receptor: Others Boiling pt: Melting pt: 175°C Solubility: DMSO: 10 mM Appearance: Off-white Color Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Sarkar SK, et al. PLoS One. 2012;7(11):e48598.

Note: This product is for in vitro research use only