

## Etofylline

Catalog No: #S0829


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Etofylline
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	224.22
Formulation	C9H12N4O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Respiratory systemSMILES:CN1C2=C(N(CCO)C=N2)C(=O)N(C)C1=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:163CSolubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Nirogi RV, et al. J Chromatogr B Analyt Technol Biomed Life Sci. 2007 Apr 1;848(2):271-6.

Note: This product is for in vitro research use only