

Benzoic acid


Catalog No: #S0833

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Benzoic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	122.12
Formulation	C7H6O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OC(=O)C1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: 249°C
Melting pt: 121-125°C
Solubility: Ethanol: 24 mg/mL (196.52 mM); DMSO: 24 mg/mL (196.52 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ohtsuki S, et al. J Neurochem. 2002 Oct;83(1):57-66.

Note: This product is for in vitro research use only