

Ethylparaben

Catalog No: #S0834


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ethylparaben
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	166.17
Formulation	C ₉ H ₁₀ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CCOC(=O)C1=CC=C(O)C=C1 Pathways: Others Receptor: Others Boiling pt: 297-298°C Melting pt: 114-117°C Solubility: DMSO: > 10 mM Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jianmei C, et al. Microbiol Res. 2015 Mar; 172:48-56.

Note: This product is for in vitro research use only