

## Triclabendazole

Catalog No: #S0838


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Triclabendazole
Brief Description	Inhibitors
Purification	95.00%
Target Name	Microtubule Associated inhibitor
Calculated MW	359.66
Formulation	C <sub>14</sub> H <sub>9</sub> Cl <sub>3</sub> N <sub>2</sub> O <sub>2</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: CSC1=NC2=CC(OC3=C(Cl)C(Cl)=CC=C3)=C(Cl)C=C2N1 Pathways: Cytoskeletal

Signaling Receptor: Tubulin Binding pt: 495.9 kDa Melting pt: 175-176°C Solubility: DMSO: 72 mg/mL (200.18 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Moll L, et al. Vet Parasitol, 2000 Jul 24, 91(1-2), 153-158.

Note: This product is for in vitro research use only