

CYP3A7 Polyclonal Antibody

Catalog No: #40818

Package Size: #40818-1 50ul #40818-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CYP3A7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CYP3A7 Polyclonal Antibody detects endogenous levels of CYP3A7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CYP3A7.
Target Name	CYP3A7
Other Names	CYP3A7; Cytochrome P450 3A7; CYP11A7; Cytochrome P450-HFLA
Accession No.	Swiss-Prot: P24462NCBI Gene ID: 1551
Uniprot	P24462
GeneID	100861540;1551;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using CYP3A7 Polyclonal Antibody

Note: This product is for in vitro research use only