

## Ademine(Triamterene)

Catalog No: #S0843


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ademine(Triamterene)
Brief Description	Inhibitors
Purification	99.00%
Target Name	Sodium Channel inhibitor
Calculated MW	253.26
Formulation	C12H11N7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular SystemSMILES:NC1=NC2=C(N=C(C(N)=N2)C2=CC=CC=C2)C(N)=N1Pathways:Membrane transporter/Ion channelReceptor:ENaC; ENaCaS583Boiling pt:Melting pt:316CSolubility:DMSO: 20 mg/mL (78.97 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Busch AE, et al. Pflugers Arch, 1996, 432(5), 760-766.

Note: This product is for in vitro research use only