

Levodopa


Catalog No: #S0848

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Levodopa
Brief Description	Inhibitors
Purification	99.64%
Target Name	Dopamine Receptor agonist
Calculated MW	197.19
Formulation	C9H11NO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: C1=CC(=C(C=C1C[C@@H](C(=O)O)N)O)O
Pathways: GPCR/G Protein;
Neuroscience
Receptor: Dopamine
Boiling pt: Melting pt: 276-278°C
Solubility: Ethanol:

References

1. Onofri M, et al. Expert Opin Investig Drugs. 2008 Jul; 17(7):1115-25.

Note: This product is for in vitro research use only