

## Pregnenolone

Catalog No: #S0851


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pregnenolone
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR agonist
Calculated MW	316.48
Formulation	C21H32O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: CC(=O)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC=C4[C@@]3(CC[C@@H](C4)O)C)C

Pathways: Neuroscience Receptor: mAChR Boiling pt: Melting pt: 188-190°C Solubility: Ethanol: 22 mg/mL (69.51 mM); DMSO: 22 mg/mL (69.51

mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wong P, et al. Transl Psychiatry. 2015 Mar 17; 5:e528.

Note: This product is for in vitro research use only