

Carbazochrome sodium sulfonate

Catalog No: #S0852


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carbazochrome sodium sulfonate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	376.32
Formulation	C10H17N4NaO8S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

SystemSMILES:CN1C(CC2=C/C(=N)NC(=O)N)/C(=O)C=C21)S(=O)(=O)[O-].[Na+]Pathways:OthersReceptor:OthersBoling pt:Melting

pt:Solubility:Water: 3 mg/mL (9.3 mM); DMSO: 64 mg/mL (198.59 mM)Appearance:Orange PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sendo T, et al. Naunyn Schmiedebergs Arch Pharmacol, 2003, 368(3), 175-180.

Note: This product is for in vitro research use only