

(R)-Naproxen


Catalog No: #S0855

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	(R)-Naproxen
Brief Description	Inhibitors
Purification	99.91%
Target Name	COX inhibitor
Calculated MW	230.26
Formulation	C ₁₄ H ₁₄ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system SMILES: C[C@H](C1=CC2=C(C=C1)C=C(C=C2)OC)C(=O)O Pathways: Immunology/Inflammation;

Neuroscience Receptor: COX Boiling pt: 403.9°C Melting pt: 156-158°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Albertini R, et al. *Inflamm Res.* 2007, 56(6):228-9.

Note: This product is for in vitro research use only