

Dehydroepiandrosterone acetate

Catalog No: #S0856


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dehydroepiandrosterone acetate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Androgen Receptor agonist
Calculated MW	330.46
Formulation	C21H30O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC(=O)O[C@H]1CC[C@@]2([C@H]3CC[C@]4([C@H]([C@@H]3CC=C2C1)CCC4=O)C)CPathways:Endocrinology/HormonesReceptor:Androgen ReceptorBoiling pt:Melting pt:168-170CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Strka L, et al. J Steroid Biochem Mol Biol. 2015 Jan;145:254-60.

Note: This product is for in vitro research use only