

Fenipentol


Catalog No: #S0857

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fenipentol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Lipid inhibitor
Calculated MW	164.25
Formulation	C ₁₁ H ₁₆ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: CCCCC(O)C1=CC=CC=C1
Pathways: Metabolism
Receptor: CES1
Boiling pt: 245.2°C
Melting pt: -20.5°C
Solubility: DMSO: 10 mM
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Chey WY, et al. *Gastroenterology*. 1983 Jun;84(6):1578-84.

Note: This product is for in vitro research use only