

Dendrid

Catalog No: #S0863


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dendrid
Brief Description	Inhibitors
Purification	99.45%
Target Name	Antifection inhibitor
Calculated MW	354.1
Formulation	C ₉ H ₁₁ N ₂ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: C1[C@@H]([C@H](O[C@H]1N2C=C(C(=O)NC2=O)I)CO)O Pathways: Microbiology&virology Receptor: Feline herpesvirus type-1 (FHV-1) Boiling pt: Melting pt: 194°C Solubility: DMSO: 35 mg/mL (98.84 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Silvestri DL, et al. JAMA. 1982 Aug 27; 248(8):953-9.

Note: This product is for in vitro research use only