

Noctone


Catalog No: #S0865

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Noctone
Brief Description	Inhibitors
Purification	99.13%
Target Name	Histamine Receptor antagonist
Calculated MW	350.86
Formulation	C13H23CIN4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: CN/C(=C/[N+](=O)[O-])/NCCSCC1=CC=C(O1)CN(C)C
 Pathways: GPCR/G Protein;
 Immunology/Inflammation; Neuroscience
 Receptor: H2 receptor
 Boiling pt: Melting pt: 134°C
 Solubility: Water: 70 mg/mL (199.5 mM); DMSO: 70 mg/mL (199.5 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Pattichis K, Louca LL. Drug Metabol Drug Interact. 1995; 12(1):1-36.

Note: This product is for in vitro research use only