

Bufexamac

Catalog No: #S0867


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bufexamac
Brief Description	Inhibitors
Purification	95.99%
Target Name	COX inhibitor
Calculated MW	223.27
Formulation	C12H17NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CCCCOC1=CC=C(CC(=O)NO)C=C1Pathways:Immunology/Inflammation;

NeuroscienceReceptor:IFN- α Boling pt:Melting pt:161-162CSolubility:Ethanol: 3 mg/mL (13.43 mM); DMSO: 45 mg/mL (201.54 mM)Appearance:White schistose crystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kr?nke B, et al. Contact Dermatitis. 1997 Apr;36(4):212-5.

Note: This product is for in vitro research use only