

## Altrenogest

Catalog No: #S0870


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Altrenogest
Brief Description	Inhibitors
Purification	99.50%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	310.43
Formulation	C21H26O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:C[C@]12C=CC3=C4CCC(=O)C=C4CC[C@H]3[C@@H]1CC[C@]2(CC=C)OPathways:Endocrinology/HormonesReceptor:ERBoiling pt:495.6C?Melting pt:120CSolubility:Ethanol: 2 mg/mL (6.44 mM); DMSO: 62 mg/mL (199.72 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Squires EL, et al. J Anim Sci, 1983, 56(4), 901-910.

Note: This product is for in vitro research use only