

Noradrenaline bitartrate

Catalog No: #S0871


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Noradrenaline bitartrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adrenergic Receptor
Calculated MW	319.26
Formulation	C8H11NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
 SMILES: OC(C(O)C(O)=O)C(O)=O.NCC(O)C1=CC(O)=C(O)C=C1
 Pathways: Neuroscience; GPCR/G Protein Receptor: Adrenergic Receptor
 Boling pt: Melting pt: 100-106°C
 Solubility: Soluble in Water
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lachnit WG, et al. Br J Pharmacol. 1997 Mar;120(5):819-26.

Note: This product is for in vitro research use only