

Amprolium

Catalog No: #S0872


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amprolium
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	278.79
Formulation	C ₁₄ H ₁₉ N ₄ Cl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: [Cl-].CCCC1=NC=C(C[N+]2=CC=CC=C2)C(N)=N1
 Pathways: Others
 Receptor: thiamine
 Boiling pt: Melting pt: Solubility: Water: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Iqbal A, et al. J Parasit Dis. 2013 Apr;37(1):88-93.

Note: This product is for in vitro research use only