

Amprolium Hydrochloride

Catalog No: #S0872L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amprolium Hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	315.24
Formulation	C ₁₄ H ₁₉ ClN ₄ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: Cl.[Cl-].CCCC1=NC=C(C[N+]2=CC=CC=C2C)C(N)=N1 Pathways: Others Receptor: thiamine Boiling pt: Melting pt: 248-249°C Solubility: Water: 63 mg/mL (199.84 mM); DMSO: 2 mg/mL (6.34 mM) Appearance: White or Off-white Color Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hamamoto K, et al. J Vet Pharmacol Ther, 2000, 23(1), 9-14.

Note: This product is for in vitro research use only