

Dimetridazole

Catalog No: #S0874


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dimetridazole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	141.13
Formulation	C5H7N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CN1C(C)=NC=C1[N+](=O)[O-]
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 177-182°C
Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Dalton P. Vet Rec. 1996 Oct 19;139(16):399.

Note: This product is for in vitro research use only