

Clioquinol

Catalog No: #S0876


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Clioquinol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	305.5
Formulation	C ₉ H ₅ ClINO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: OC1=C(I)C=C(Cl)C2=C1N=CC=C2
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 175-183°C
Solubility: Chloroform: 15mg/mL
Appearance: Faint Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Rohde W, et al. Antimicrob Agents Chem Others. 1976 Aug;10(2):234-40.

Note: This product is for in vitro research use only