

## Gallic acid

Catalog No: #S0877


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Gallic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	170.12
Formulation	C7H6O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: OC(=O)C1=CC(O)=C(O)C(O)=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: 252°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Liu Y, et al. Bioorg Med Chem Lett. 2013 Dec 1;23(23):6336-40.

Note: This product is for in vitro research use only