

Diosmin

Catalog No: #S0878


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diosmin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Norepinephrine
Calculated MW	608.54
Formulation	C28H32O15
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)OC[C@@H]2[C@H]([C@@H]([C@H]([C@@H](O2)OC3=CC(=C4C(=C3)OC(=CC4=O)C5=CC(=C(C=C5)OC)O)O)O)O)O)OPathways:NeuroscienceReceptor:Norepinephrine (NE)Boiling pt:Melting pt:277-278CSolubility:DMSO: 6 mg/mL (9.85 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Veverkov L, et al. Rozhl Chir. 2005 Aug;84(8):410-2, 414-6.

Note: This product is for in vitro research use only